# Spectral\_image class guide

isabelpostmes

January 2021

#### 1 Guide

#### 1.1 Loading data

Lets talk through the Spectral\_image class. We start by loading a spectral image, saved in a .dm3 or .dm4 file through:

```
>>> im = Spectral_image.load_data('path/to/dmfile.dm4') 1
```

This calls on an alternative constructor, in which the data from the dm-file is loaded, and plugged into the regular constructor. In this function, the loading package ncempy.io.dm is used, more info here.

```
81
        @classmethod
82
        def load_data(cls, path_to_dmfile):
83
            INPUT:
84
                path_to_dmfile: str, path to spectral image file (.
85
                    dm3 or .dm4 extension)
86
            OUTPUT:
                image — Spectral_image, object of Spectral_image
87
                    class containing the data of the dm-file
88
            dmfile = dm.fileDM(path_to_dmfile).getDataset(0)
89
            data = np.swapaxes(np.swapaxes(dmfile['data'], 0,1), 1,
90
                2)
91
            ddeltaE = dmfile['pixelSize'][0]
92
            pixelsize = np.array(dmfile['pixelSize'][1:])
93
            energyUnit = dmfile['pixelUnit'][0]
94
            ddeltaE *= cls.get_prefix(energyUnit, 'eV')
            pixelUnit = dmfile['pixelUnit'][1]
95
            pixelsize *= cls.get_prefix(pixelUnit, 'm')
97
            image = cls(data, ddeltaE, pixelsize = pixelsize)
            return image
98
```

Furthermore, we see the cls.get\_prefix(), which is a small function which recognises the prefix in a unit and transfers it to a numerical value (e.g. 1E3 for k), see lines 870-916 in the complete code. Furthermore, the general constructor is called upon with cls(data, ddeltaE, pixelsize = pixelsize).

The spectral image class starts by defining some constant variables, both class related and physical, which you can find in the complete code. The class constructor takes in at least the data of the spectral image, data, and the broadness of the energy loss bins, deltadeltaE. Other metadata can be given if known.

Also, the delta\_E axis, that is the energy-loss axis is determined by self.determine\_deltaE(), based upon the broadness of the energy-loss bins and the index at which the average of all spectra in the image has it maximum. The definition of determine\_deltaE() can be found at line 101 in the complete code. The image axes are determined by self.calc\_axes(), and output either index arrays, or, if the pixel size is defined, the spacings array in meters. The definition of calc\_axes() can be found at line 116 in the complete code.

The definitions of some other properties, such as im.1, im.image\_shape, and im.shape can be found in the complete code from line 69 onwards.

Furthermore, there are some retrieving functions, such as im.get\_data(), im.get\_deltaE(), im.get\_metadata, and get\_pixel\_signal, which should be quite self-explainatory, but whose definitions can be found in the complete code from line 124 onwards.

```
53
        def __init__(self, data, deltadeltaE, pixelsize = None,
           beam_energy = None, collection_angle = None, name = None
           ):
54
            self.data = data
55
            self.ddeltaE = deltadeltaE
56
            self.determine_deltaE()
57
            if pixelsize is not None:
58
                self.pixelsize = pixelsize
59
            self.calc_axes()
```

self.beam\_energy = beam\_energy

if collection\_angle is not None:

if beam\_energy is not None:

if name is not None:

self.name = name

# 1.2 Preparing data

60 61

62

63

64

65

class Spectral\_image():

Now that we have loaded the data, we can perform some operations on the data of the image before starting any calculations, if wished. For example, if we wish to cut the image to a rectangle ranging from pixel a trough pixel b in

self.collection\_angle = collection\_angle

width and from pixel c trough pixel d in height, you can simply run (added the +1's to emphasise the excluding nature of the function, for definition of function cut\_image, see line 154):

```
>>> im.cut_image([a,b+1], [c,d+1])
```

Also, one can, in the future, cut to a certain energy range from E1 to E2, by running (for definition of function cut, see line 150):

```
>>> im.cut(E1, E2)
```

Also, one can decide to smooth the spectra. The default smoothing is done by convoluting a length 10 Hanning window, but this can be altered by adding arguments window\_len= , and window= respectively to the call function. Also it should be noted that by defualt, the original spectra are disregarded and overwritten by the smoothed signal, to save memory. If you do not want this, add keep\_original=True to your call-statement. Please note that in this case to call upon your smoothed data, you should call im.data\_smooth in stead of simply im.data. The definition of function smooth can be found in the complete code from line 164. For smoothing your data by a moving average over 50 values for example, run:

```
>>> im.smooth(window_len= 50, window = 'flat')
```

### 1.3 Calculations on image

Now that you have altered your image to your wishes, we can start the calculations on the spectra.

One of the first things you probably want to do (otherwise why are you using this class instead of hyperspy), is loading in the trained ZLPs for the image. This can be done by running (for the definition of function calc\_ZLPs\_gen2, see complete code from line 446):

```
>>> im.calc_ZLPs_gen2()
```

Due to memory considerations, it is not adviced to calculate the ZLPs for each pixel at once, but calculate them per pixel as they are needed. This can be done for pixel at coordinate [i,j] by (for the definition of function calc\_ZLPs, see complete code from line 446):

```
>>> ZLPs_pixel_i_j = im.calc_ZLPs(i,j)
```

NB: the two functions mentioned above are at the time mere near copies from the code of Laurien, they might/will change significantly when the training of the neural network for ZLPs is incorperated into this class.

When the ZLPs are calculated, one can deconvolute the EEL spectra, to obtain the single scattering distributions. This is done as Egerton explains in his book, and the definition of function calc\_ZLPs can be found in the complete code from line 202 on. This is also done per pixel and per ZLP, once again for memory considerations.

```
>>> S_E_ijk = im.deconvolute(i,j,ZLP_k)
```

Both the calculation of the ZLPs and the deconvolution of the signal are needed for the evaluation with the Kramers Kronig analysis. The Kramers Kronigs analysis inplemented in our code is the version Egerton explained in his book, with the hyperspy adaptation for tail correction. The definition of function kramers\_kronig\_hs can be found in the complete code from line 481 on. For a single scattering distribution S\_E\_ijk and integrated ZLP\_k intensity N\_ZLP\_k, you can call get the dielectric function, thickness, and an approximation of the surface scattering contribution in the single scattering distribution by running:

```
>>> df, t, SS_E = im.kramers_kronig_hs(S_E_ijk, N_ZLP_k)
```

If one wishes to calculate for each pixel the average dielectric function and thickness, you can call upon the im\_dielectric\_function function, whose definition can be found at line 681. It subsequently calculates for each pixel [i,j], for each predicited ZLP\_k, the dielectric function and thickness, and saves as an attribute for each pixel the average dielectric function and average thickness, and the standard deviation in both.

```
>>> im.im_dielectric_function()
```

Now that the dielectric function is calculated, we can evaluate it, for example by considering the crossings of the real part of the dielectric function from negative to positive. If the function crossings\_im (line 756) is called upon, there are two attributes created: im.crossings\_E and im.crossings\_n, where in the first all the energy values at which the average dielectric function of each pixel crosses are saved, and in the latter the number of crossings at each pixel.

```
>>> im.crossings_im()
```

# 1.4 Plotting functions

At this point, two plotting functions are implemented: plot\_sum (line 804) and plot\_all (line 835). The first plots the integrated intensity at each pixel.

```
>>> im.plot_sum()
```

With the latter, you can plot all spectra in the image, either in a single plot (default), or in a single image each (set same\_image =False). Furthermore you can choose a range of pixels for the spectra you want to plot, the range of energy, and change from the IEELS default, to plotting other functions per pixel, such as the average dielectric function. Other functionels can be found in the definition (line 835 complete code).

```
>>> im.plot_all(self, range_x = [10,20], range_y = [0,80],
    range_E = [0.8,5], signal = "dielectric_function")
```

# 2 Complete code

```
1 #!/usr/bin/env python3
2 # -*- coding: utf-8 -*-
3
4 Created on Fri Dec 18 17:38:57 2020
5
6 @author: isabel
7
8 #!/usr/bin/env python3
9 # -*- coding: utf-8 -*-
  0.00
10
11 Created on Thu Dec 17 23:12:19 2020
12
13 @author: isabel
14 """
15 #!/usr/bin/env python3
16 # -*- coding: utf-8 -*-
17 """
18 Created on Tue Nov 10 01:05:56 2020
19
20 @author: isabel
21
22
23 import pandas as pd
24 import glob
25 import matplotlib.pyplot as plt
26 import seaborn as sns
27 import natsort
28 import tensorflow.compat.v1 as tf
29 import seaborn as sns
30 import numpy as np
31 #from lmfit import Model
32 from scipy.fftpack import next_fast_len
33 import logging
34 from ncempy.io import dm;
35
36 tf.get_logger().setLevel('ERROR')
37
39
41 class Spectral_image():
       DIELECTRIC_FUNCTION_NAMES = ['dielectric_function', 'dielectricfunction', '
           dielec_func', 'die_fun', 'df', 'epsilon']
```

```
43
        EELS_NAMES = ['electron_energy_loss_spectrum','electron_energy_loss','EELS', 'EEL',
             'energy_loss', 'data']
44
        IEELS_NAMES = ['inelastic_scattering_energy_loss_spectrum', '
            inelastic_scattering_energy_loss', 'inelastic_scattering', 'IEELS', 'IES']
45
        ZLP_NAMES = ['zeros_loss_peak', 'zero_loss', 'ZLP', 'ZLPs']
46
        m_{-}0 = 511.06 \text{ #eV}, electron rest mass
47
        a_0 = 5.29E-11 \text{ #m}, Bohr radius
48
        h_bar = 6.582119569E-16 \#eV/s
49
50
        c = 2.99792458E8 \#m/s
51
52
53
        def __init__(self, data, deltadeltaE, pixelsize = None, beam_energy = None,
            collection_angle = None, name = None):
54
            self.data = data
55
            self.ddeltaE = deltadeltaE
56
            self.determine_deltaE()
57
            if pixelsize is not None:
58
                self.pixelsize = pixelsize
59
            self.calc_axes()
60
            if beam_energy is not None:
61
                self.beam_energy = beam_energy
62
            if collection_angle is not None:
63
                self.collection_angle = collection_angle
64
            if name is not None:
65
                self.name = name
66
67
        #PROPERTIES
68
69
        @property
70
        def l(self):
71
            return self.data.shape[2]
72
        @property
73
        def image_shape(self):
74
            return self.data.shape[:2]
75
76
        @property
77
        def shape(self):
78
            self.shape = self.data.shape
79
80
81
        @classmethod
        def load_data(cls, path_to_dmfile):
82
83
            INPUT:
84
                path_to_dmfile: str, path to spectral image file (.dm3 or .dm4 extension)
85
```

```
86
             OUTPUT:
87
                 image -- Spectral_image, object of Spectral_image class containing the data
                      of the dm-file
88
89
             dmfile = dm.fileDM(path_to_dmfile).getDataset(0)
90
             data = np.swapaxes(np.swapaxes(dmfile['data'], 0,1), 1,2)
91
             ddeltaE = dmfile['pixelSize'][0]
             pixelsize = np.array(dmfile['pixelSize'][1:])
92
93
             energyUnit = dmfile['pixelUnit'][0]
             ddeltaE *= cls.get_prefix(energyUnit, 'eV')
94
             pixelUnit = dmfile['pixelUnit'][1]
95
96
             pixelsize *= cls.get_prefix(pixelUnit, 'm')
97
             image = cls(data, ddeltaE, pixelsize = pixelsize)
98
             return image
99
100
101
        def determine_deltaE(self):
102
             INPUT:
103
104
                 self
105
             Determines the delta energies of the spectral image, based on the delta delta
106
                 energie,
107
             and the index on which the spectral image has on average the highest intesity,
                this
108
             is taken as the zero point for the delta energy.
109
110
             data_avg = np.average(self.data, axis = (0,1))
111
             ind_max = np.argmax(data_avg)
112
             self.deltaE = np.linspace(-ind_max * self.ddeltaE, (self.l-ind_max-1)*self.
                 ddeltaE, self.1)
             #return deltaE
113
114
115
116
        def calc_axes(self):
117
             self.y_axis = np.linspace(0, self.image_shape[0]-1, self.image_shape[0])
             self.x_axis = np.linspace(0, self.image_shape[1]-1, self.image_shape[1])
118
             if hasattr(self, 'pixelsize'):
119
120
                 self.y_axis *= self.pixelsize[0]
121
                 self.x_axis *= self.pixelsize[1]
122
        #RETRIEVING FUNCTIONS
123
124
        def get_data(self):
125
             return self.data
126
        def get_deltaE(self):
127
```

```
128
             return self.deltaE
129
130
        def get_metadata(self):
131
             meta_data = \{\}
132
             if self.beam_energy is not None:
133
                 meta_data['beam_energy'] = self.beam_energy
134
             if self.collection_angle is not None:
135
                 meta_data['collection_angl'] = self.collection_angle
136
             return meta_data
137
138
        def get_pixel_signal(self, i,j, signal = 'EELS'):
             #TODO: add alternative signals + names
139
140
             if signal == 'EELS':
141
                 return np.copy(self.data[ i, j, :])
142
             elif signal == 'df_avg':
143
                 return np.copy(self.dielectric_function_im_avg[ i, j, :])
144
             else:
145
                 return np.copy(self.data[ i, j, :])
146
147
        #METHODS ON SIGNAL
148
149
150
        def cut(self, E1, E2):
             #TODO
151
152
             pass
153
        def cut_image(self, range_width, range_height):
154
155
             #TODO: add floats for cutting to meter sizes?
             self.data = self.data[range_height[0]:range_height[1], range_width[0]:
156
                 range_width[1]]
157
             self.y_axis = self.y_axis[range_height[0]:range_height[1]]
158
             self.x_axis = self.x_axis[range_width[0]:range_width[1]]
159
160
        #TODO
161
        def samenvoegen(self):
162
             pass
163
164
        def smooth(self, window_len=10, window='hanning', keep_original = False):
165
             """smooth the data using a window with requested size.
166
             This method is based on the convolution of a scaled window with the signal.
167
168
             The signal is prepared by introducing reflected copies of the signal
             (with the window size) in both ends so that transient parts are minimized
169
             in the begining and end part of the output signal.
170
171
172
             input:
```

```
173
                 x: the input signal
174
                 window_len: the dimension of the smoothing window; should be an odd integer
175
                 window: the type of window from 'flat', 'hanning', 'hamming', 'bartlett', '
                     blackman'
176
                     flat window will produce a moving average smoothing.
177
178
             output:
179
                 the smoothed signal
180
             .....
181
             #TODO: add comnparison
182
             window_len += (window_len+1)%2
183
184
             s=np.r_{['-1', self.data[:,:,window_len-1:0:-1], self.data, self.data[:,:,-2:-1]}
                 window_len-1:-1]
185
186
             if window == 'flat': #moving average
187
                 w=np.ones(window_len,'d')
188
             else:
189
                 w=eval('np.'+window+'(window_len)')
190
191
             #y=np.convolve(w/w.sum(),s,mode='valid')
             surplus_data = int((window_len-1)*0.5)
192
193
             if keep_original:
194
                 self.data_smooth = np.apply_along_axis(lambda m: np.convolve(m, w/w.sum(),
                     mode='valid'), axis=2, arr=s)[:,:,surplus_data:-surplus_data]
195
             else:
196
                 self.data = np.apply_along_axis(lambda m: np.convolve(m, w/w.sum(), mode='
                     valid'), axis=2, arr=s)[:,:,surplus_data:-surplus_data]
197
198
199
             return #y[(window_len-1):-(window_len)]
200
201
202
        def deconvolute(self, i,j, ZLP):
203
204
             y = self.get_pixel_signal(i,j)
205
             r = 3 #Drude model, can also use estimation from exp. data
206
             A = y[-1]
207
             n_{times_{extra}} = 2
208
             sem_inf = next_fast_len(n_times_extra*self.1)
209
210
211
212
             y_extrp = np.zeros(sem_inf)
213
             y_ZLP_extrp = np.zeros(sem_inf)
```

```
214
             x_{extrp} = np.linspace(self.deltaE[0] - self.l*self.ddeltaE, sem_inf*self.ddeltaE
                 +self.deltaE[0] - self.l*self.ddeltaE, sem_inf)
215
             x_{extrp} = np.linspace(self.deltaE[0], sem_inf*self.ddeltaE+self.deltaE[0],
216
                 sem_inf)
217
             y_ZLP_extrp[:self.1] = ZLP
218
219
             y_extrp[:self.1] = y
             x_extrp[:self.1] = self.deltaE[-self.1:]
220
221
             y_extrp[self.l:] = A*np.power(1+x_extrp[self.l:]-x_extrp[self.l],-r)
222
223
224
             x = x_extrp
225
             y = y_extrp
226
             y_ZLP = y_ZLP_extrp
227
228
             z_nu = CFT(x,y_ZLP)
229
             i_nu = CFT(x,y)
230
             abs_i_nu = np.absolute(i_nu)
231
             N_{ZLP} = 1 + scipy.integrate.cumtrapz(y_{ZLP}, x, initial = 0)[-1] + 1 + arbitrary units
                 ??? np.sum(EELZLP)
232
233
             s_nu = N_ZLP*np.log(i_nu/z_nu)
234
             j_{-nu} = z_{-nu} * s_{-nu} / N_{-}ZLP
235
             S_E = np.real(iCFT(x,s_nu))
236
             s_nu_nc = s_nu
             s_nu_nc[500:-500] = 0
237
238
             S_E_nc = np.real(iCFT( x,s_nu_nc))
239
             J_{L}^{T} = np.real(iCFT(x,j_{nu}))
240
241
             return J1_E[:self.1]
242
243
         #METHODS ON ZLP
244
         #CALCULATING ZLPs FROM PRETRAINDED MODELS
245
         def calculate_general_ZLPs(self, path_to_models):
246
             tf.reset_default_graph()
247
             #TODO: redifine paths based upon new fitter saving modes
             #TODO: rewrite to have models as atributes?
248
249
250
             d_{string} = '07.09.2020'
             path_to_data = 'Data_oud/Results/%(date)s/'% {"date": d_string}
251
252
             path_predict = r'Predictions_*.csv'
253
254
             path_cost = r'Cost_*.csv'
255
256
             all_files = glob.glob(path_to_data + path_predict)
```

```
257
258
             li = []
259
             for filename in all_files:
260
                 df = pd.read_csv(filename, delimiter=",", header=0, usecols=[0,1,2], names
                     =['x', 'y', 'pred'])
261
                 li.append(df)
262
263
             training_data = pd.concat(li, axis=0, ignore_index=True)
264
265
             self.dE1 = np.round(max(training_data['x'][(training_data['x']< 3)]),2)</pre>
266
             self.dE2 = np.round(min(training_data['x'][(training_data['x']> 3)]),1)
267
             self.dE0 = np.round(self.dE1 - .5, 2)
268
269
             all_files_cost = glob.glob(path_to_data + path_cost)
270
             all_files_cost_sorted = natsort.natsorted(all_files_cost)
271
272
             chi2_array = []
273
             chi2_index = []
274
275
             for filename in all_files_cost_sorted:
276
                 df = pd.read_csv(filename, delimiter=",", header=0, usecols=[0,1], names=['
                     train', 'test'])
277
                 best_try = np.argmin(df['test'])
278
                 chi2_array.append(df.iloc[best_try,0])
279
                 chi2_index.append(best_try)
280
281
             chi_data = pd.DataFrame()
282
             chi_data['Best chi2 value'] = chi2_array
             chi_data['Epoch'] = chi2_index
283
284
285
             good_files = []
286
             count = 0
287
             threshold = 3
288
289
             for i,j in enumerate(chi2_array):
290
                 if j < threshold:</pre>
291
                     good_files.append(1)
292
                     count += 1
293
                 else:
294
                     good_files.append(0)
295
296
             tf.get_default_graph()
297
             tf.disable_eager_execution()
             #config = tf.ConfigProto()
298
299
             #config.gpu_options.allow_growth = True
300
```

```
301
302
             def make_model(inputs, n_outputs):
                 hidden_layer_1 = tf.layers.dense(inputs, 10, activation=tf.nn.sigmoid)
303
304
                 hidden_layer_2 = tf.layers.dense(hidden_layer_1, 15, activation=tf.nn.
                     sigmoid)
                hidden_layer_3 = tf.layers.dense(hidden_layer_2, 5, activation=tf.nn.relu)
305
306
                 output = tf.layers.dense(hidden_layer_3, n_outputs, name='outputs', reuse=
                     tf.AUTO_REUSE)
307
                 return output
308
309
            x = tf.placeholder("float", [None, 1], name="x")
310
             predictions = make_model(x, 1)
311
312
313
             prediction_file = pd.DataFrame()
314
             len_data = self.1
315
             predict_x = np.linspace(-0.5, 20, 1000).reshape(1000,1)
316
             predict_x = self.deltaE.reshape(self.1,1)
317
318
             self.ZLPs_gen = np.zeros((count, len_data))
319
             with tf.Session() as sess: #TODO: gives warning
320
                 sess.run(tf.global_variables_initializer())
321
322
                 for i in range(0,len(good_files)):
323
                     if good_files[i] == 1:
324
                         best_model = 'Models_oud/Best_models/%(s)s/best_model_%(i)s'% {'s':
                              d_string, 'i': i}
325
                         saver = tf.train.Saver(max_to_keep=1000)
326
                         saver.restore(sess, best_model)
327
328
                         extrapolation = sess.run(predictions, #TODO: RESTARTS KERNEL!!!!!
329
                                                  feed_dict={
330
                                                  x: predict_x
331
                                                  })
332
                         prediction_file['prediction_%(i)s' % {"i": i}] = extrapolation.
                             reshape(len_data,)
333
                         self.ZLPs_gen[i, :] = np.exp(extrapolation)#.reshape(len_data,)
334
335
336
        @staticmethod
337
        def make_model(inputs, n_outputs):
338
             hidden_layer_1 = tf.layers.dense(inputs, 10, activation=tf.nn.sigmoid)
339
             hidden_layer_2 = tf.layers.dense(hidden_layer_1, 15, activation=tf.nn.sigmoid)
             hidden_layer_3 = tf.layers.dense(hidden_layer_2, 5, activation=tf.nn.relu)
340
341
             output = tf.layers.dense(hidden_layer_3, n_outputs, name='outputs', reuse=tf.
                AUTO_REUSE)
```

```
342
             return output
343
344
        def calc_ZLPs_gen2(self, specimen = 4):
345
             tf.reset_default_graph()
346
             if specimen == 3:
                 d_{string} = '06.12.2020'
347
                 path_to_data = 'Data_oud/Results/sp3/%(date)s/'% {"date": d_string}
348
349
             else:
350
                 d_string = '07.09.2020'
351
                 path_to_data = 'Data_oud/Results/%(date)s/'% {"date": d_string}
352
353
             path_predict = r'Predictions_*.csv'
354
             path_cost = r'Cost_*.csv'
355
             all_files = glob.glob(path_to_data + path_predict)
356
357
358
             li = []
359
             for filename in all_files:
                 df = pd.read_csv(filename, delimiter=",", header=0, usecols=[0,1,2], names
360
                     =['x', 'y', 'pred'])
361
                 li.append(df)
362
363
             training_data = pd.concat(li, axis=0, ignore_index=True)
364
365
366
367
             all_files_cost = glob.glob(path_to_data + path_cost)
368
369
370
             import natsort
371
372
             all_files_cost_sorted = natsort.natsorted(all_files_cost)
373
374
             chi2_array = []
375
             chi2_index = []
376
377
             for filename in all_files_cost_sorted:
                 df = pd.read_csv(filename, delimiter=",", header=0, usecols=[0,1], names=['
378
                     train', 'test'])
379
                 best_try = np.argmin(df['test'])
380
                 chi2_array.append(df.iloc[best_try,0])
381
                 chi2_index.append(best_try)
382
383
             chi_data = pd.DataFrame()
             chi_data['Best chi2 value'] = chi2_array
384
385
             chi_data['Epoch'] = chi2_index
```

```
386
387
388
389
             good_files = []
390
             count = 0
             threshold = 3
391
392
393
             for i,j in enumerate(chi2_array):
                 if j < threshold:</pre>
394
395
                      good_files.append(1)
                      count += 1
396
397
                 else:
398
                      good_files.append(0)
399
400
401
402
403
             tf.get_default_graph
404
             tf.disable_eager_execution()
405
406
407
             x = tf.placeholder("float", [None, 1], name="x")
408
             predictions = self.make_model(x, 1)
409
410
411
             prediction_file = pd.DataFrame()
412
             len_data = self.1
413
             predict_x = np.linspace(-0.5, 20, 1000).reshape(1000, 1)
414
415
             predict_x = self.deltaE.reshape(len_data,1)
416
417
             self.ZLPs_gen = np.zeros((count, len_data))
418
             j = 0
419
             with tf.Session() as sess:
420
                 sess.run(tf.global_variables_initializer())
421
422
                 for i in range(0,len(good_files)):
423
                      if good_files[i] == 1:
424
                          if specimen ==3:
425
                              best_model = 'Models_oud/Best_models/sp3/%(s)s/best_model_%(i)s
                                  '% {'s': d_string, 'i': i}
426
                          else:
427
                              best_model = 'Models_oud/Best_models/%(s)s/best_model_%(i)s'% {
                                  's': d_string, 'i': i}
428
                          saver = tf.train.Saver(max_to_keep=1000)
429
                          saver.restore(sess, best_model)
```

```
430
431
                         extrapolation = sess.run(predictions,
432
                                                  feed_dict={
433
                                                  x: predict_x
434
                                                  })
                         #prediction_file['prediction_%(i)s' % {"i": i}] = extrapolation.
435
                             reshape(1000,)
436
                         self.ZLPs_gen[j,:] = np.exp(extrapolation.reshape(len_data,))
                         prediction_file['prediction_%(i)s' % {"i": i}] = extrapolation.
437
                             reshape(len_data,)
438
                         j += 1
439
440
             self.dE1 = np.round(max(training_data['x'][(training_data['x']< 3)]),2)</pre>
441
             self.dE2 = np.round(min(training_data['x'][(training_data['x']> 3)]),1)
442
             self.dE0 = np.round(self.dE1 - .5, 2)
443
             #return ZLPs_gen, dE0, dE1, dE2
444
445
446
        def calc_ZLPs(self, i,j):
447
             ### Definition for the matching procedure
448
             signal = self.get_pixel_signal(i,j)
449
450
             if not hasattr(self, 'ZLPs_gen'):
451
                 self.calc_ZLPs_gen2("iets")
452
453
             def matching( signal, ind_ZLP):
                 gen_i_ZLP = self.ZLPs_gen[ind_ZLP, :]*np.max(signal)/np.max(self.ZLPs_gen[
454
                     ind_ZLP,:]) #TODO!!!!, normalize?
455
                 delta = np.divide((self.dE1 - self.dE0), 3)
456
457
                 factor_NN = np.exp(- np.divide((self.deltaE[(self.deltaE<self.dE1) & (self.</pre>
                     deltaE >= self.dE0)] - self.dE1)**2, delta**2))
458
                 factor_dm = 1 - factor_NN
459
460
                 range_0 = signal[self.deltaE < self.dE0]
461
                 range_1 = gen_i_ZLP[(self.deltaE < self.dE1) & (self.deltaE >= self.dE0)] *
                      factor_NN + signal[(self.deltaE < self.dE1) & (self.deltaE >= self.dE0
                     )] * factor_dm
462
                 range_2 = gen_i_ZLP[(self.deltaE >= self.dE1) & (self.deltaE < 3 * self.dE2</pre>
                     )]
463
                 range_3 = gen_i_ZLP[(self.deltaE >= 3 * self.dE2)] * 0
464
                 totalfile = np.concatenate((range_0, range_1, range_2, range_3), axis=0)
465
                 #TODO: now hardcoding no negative values!!!! CHECKKKK
                 totalfile = np.minimum(totalfile, signal)
466
467
                 return totalfile
```

468

```
469
             count = self.ZLPs_gen.shape[0]
470
             ZLPs = np.zeros(self.ZLPs_gen.shape) #np.zeros((count, len_data))
471
472
473
             for k in range(count):
                 ZLPs[k,:] = matching(signal, k)#matching(energies, np.exp(mean_k), data)
474
475
476
             return ZLPs
477
478
        #METHODS ON DIELECTRIC FUNCTIONS
479
480
481
        def kramers_kronig_hs(self, I_EELS,
482
                                 N_ZLP=None,
483
                                 iterations=1,
484
                                 n=None,
485
                                 t=None,
486
                                 delta=0.5, correct_S_s = False):
487
             r"""Calculate the complex
488
             dielectric function from a single scattering distribution (SSD) using
489
             the Kramers-Kronig relations.
490
             It uses the FFT method as in [1]. The SSD is an
491
492
             EELSSpectrum instance containing SSD low-loss EELS with no zero-loss
493
             peak. The internal loop is devised to approximately subtract the
494
             surface plasmon contribution supposing an unoxidized planar surface and
             neglecting coupling between the surfaces. This method does not account
495
496
             for retardation effects, instrumental broading and surface plasmon
497
             excitation in particles.
498
499
             Note that either refractive index or thickness are required.
500
             If both are None or if both are provided an exception is raised.
501
502
            Parameters
503
504
             zlp: {None, number, Signal1D}
                ZLP intensity. It is optional (can be None) if 't' is None and 'n'
505
                 is not None and the thickness estimation is not required. If 't'
506
                is not None, the ZLP is required to perform the normalization and
507
                 if 't' is not None, the ZLP is required to calculate the thickness.
508
509
                 If the ZLP is the same for all spectra, the integral of the ZLP
510
                 can be provided as a number. Otherwise, if the ZLP intensity is not
511
                 the same for all spectra, it can be provided as i) a Signal1D
                 of the same dimensions as the current signal containing the ZLP
512
513
                 spectra for each location ii) a BaseSignal of signal dimension 0
                 and navigation_dimension equal to the current signal containing the
514
```

```
515
                 integrated ZLP intensity.
516
             iterations: int
517
                 Number of the iterations for the internal loop to remove the
518
                 surface plasmon contribution. If 1 the surface plasmon contribution
519
                 is not estimated and subtracted (the default is 1).
520
             n: {None, float}
                 The medium refractive index. Used for normalization of the
521
522
                 SSD to obtain the energy loss function. If given the thickness
523
                 is estimated and returned. It is only required when 't' is None.
524
             t: {None, number, Signal1D}
525
                 The sample thickness in nm. Used for normalization of the
526
                  to obtain the energy loss function. It is only required when
527
                 'n' is None. If the thickness is the same for all spectra it can be
528
                 given by a number. Otherwise, it can be provided as a BaseSignal
529
                 with signal dimension 0 and navigation_dimension equal to the
530
                 current signal.
531
             delta : float
532
                 A small number (0.1-0.5 eV) added to the energy axis in
                 specific steps of the calculation the surface loss correction to
533
534
                 improve stability.
535
             full_output : bool
                 If True, return a dictionary that contains the estimated
536
                 thickness if 't' is None and the estimated surface plasmon
537
538
                 excitation and the spectrum corrected from surface plasmon
539
                 excitations if 'iterations' > 1.
540
541
             Returns
542
             eps: DielectricFunction instance
543
544
                 The complex dielectric function results,
545
546
                      .. math::
547
                          \ensuremath{\setminus} epsilon = \ensuremath{\setminus} epsilon_1 + i*\ensuremath{\setminus} epsilon_2,
548
                 contained in an DielectricFunction instance.
549
550
             output: Dictionary (optional)
551
                 A dictionary of optional outputs with the following keys:
552
553
                 "thickness"
554
                      The estimated thickness in nm calculated by normalization of
555
                      the SSD (only when 't' is None)
556
                 "surface plasmon estimation"
557
                    The estimated surface plasmon excitation (only if
558
559
                     'iterations' > 1.)
```

560

```
561
             Raises
562
563
             ValuerError
564
                 If both 'n' and 't' are undefined (None).
565
             AttribureError
566
                 If the beam_energy or the collection semi-angle are not defined in
567
                 metadata.
568
569
             Notes
570
             This method is based in Egerton's Matlab code [1] with some
571
572
             minor differences:
573
             * The wrap-around problem when computing the ffts is workarounded by
574
575
               padding the signal instead of substracting the reflected tail.
576
577
             .. [1] Ray Egerton, "Electron Energy-Loss Spectroscopy in the Electron
578
                Microscope", Springer-Verlag, 2011.
579
             ....
580
581
             output = {}
             # Constants and units
582
            me = 511.06
583
584
             e0 = 200 \# keV
585
586
             beta =30 \text{ #mrad}
587
             eaxis = self.deltaE[self.deltaE>0] #axis.axis.copy()
588
             S_E = I_EELS[self.deltaE>0]
589
590
            y = I_EELS[self.deltaE>0]
591
            1 = len(eaxis)
592
            i0 = N_ZLP
593
594
             # Kinetic definitions
             ke = e0 * (1 + e0 / 2. / me) / (1 + e0 / me) ** 2
595
596
             tgt = e0 * (2 * me + e0) / (me + e0)
597
             rk0 = 2590 * (1 + e0 / me) * np.sqrt(2 * ke / me)
598
599
             for io in range(iterations):
600
                 # Calculation of the ELF by normalization of the SSD
                 # We start by the "angular corrections"
601
602
                 Im = y / (np.log(1 + (beta * tgt / eaxis) ** 2)) / self.ddeltaE#axis.scale
603
                 if n is None and t is None:
                     raise ValueError("The thickness and the refractive index are "
604
                                       "not defined. Please provide one of them.")
605
606
                 elif n is not None and t is not None:
```

```
607
                     raise ValueError("Please provide the refractive index OR the "
608
                                       "thickness information, not both")
609
                 elif n is not None:
610
                     # normalize using the refractive index.
611
                     K = np.sum(Im/eaxis)*self.ddeltaE
612
                     K = (K / (np.pi / 2) / (1 - 1. / n ** 2))
613
                     te = (332.5 * K * ke / i0)
                 elif t is not None:
614
                     if N_ZLP is None:
615
                         raise ValueError("The ZLP must be provided when the "
616
617
                                           "thickness is used for normalization.")
                     # normalize using the thickness
618
619
                     K = t * i0 / (332.5 * ke)
620
                     te = t
621
                Im = Im / K
622
623
                 # Kramers Kronig Transform:
624
                 # We calculate KKT(Im(-1/epsilon))=1+Re(1/epsilon) with FFT
                 # Follows: D W Johnson 1975 J. Phys. A: Math. Gen. 8 490
625
626
                 # Use an optimal FFT size to speed up the calculation, and
                 # make it double the closest upper value to workaround the
627
                 # wrap-around problem.
628
                 esize = next_fast_len(2*1) #2**math.floor(math.log2(1)+1)*4
629
630
                 q = -2 * np.fft.fft(Im, esize).imag / esize
631
632
                q[:1] *= -1
                 q = np.fft.fft(q)
633
634
                 # Final touch, we have Re(1/eps)
635
                Re = q[:1].real + 1
636
                 # Egerton does this to correct the wrap-around problem, but in our
637
                 # case this is not necessary because we compute the fft on an
                 # extended and padded spectrum to avoid this problem.
638
639
                # Re=real(q)
640
                # Tail correction
641
                 # vm=Re[axis.size-1]
                 # Re[:(axis.size-1)]=Re[:(axis.size-1)]+1-(0.5*vm*((axis.size-1)) /
642
                 # (axis.size*2-arange(0,axis.size-1)))**2)
643
                 # Re[axis.size:]=1+(0.5*vm*((axis.size-1) /
644
                 # (axis.size+arange(0,axis.size)))**2)
645
646
647
                 # Epsilon appears:
648
                 # We calculate the real and imaginary parts of the CDF
649
                 e1 = Re / (Re ** 2 + Im ** 2)
650
                 e^2 = Im / (Re ** 2 + Im ** 2)
651
                if iterations > 0 and N_ZLP is not None:
652
```

```
653
                     # Surface losses correction:
654
                     # Calculates the surface ELF from a vaccumm border effect
655
                     # A simulated surface plasmon is subtracted from the ELF
656
                     Srfelf = 4 * e2 / ((e1 + 1) ** 2 + e2 ** 2) - Im
657
                     adep = (tgt / (eaxis + delta) *
658
                             np.arctan(beta * tgt / eaxis) -
659
                             beta / 1000. /
                             (beta ** 2 + eaxis ** 2. / tgt ** 2))
660
                     Srfint = 2000 * K * adep * Srfelf / rk0 / te * self.ddeltaE #axis.scale
661
662
                     if correct_S_s == True:
663
                         print("correcting S_s")
664
                         Srfint[Srfint<0] = 0</pre>
665
                         Srfint[Srfint>S_E] = S_E[Srfint>S_E]
666
                     y = S_E - Srfint
667
                     _logger.debug('Iteration number: %d / %d', io + 1, iterations)
668
669
670
             eps = (e1 + e2 * 1j)
671
             del y
672
             del I_EELS
             if 'thickness' in output:
673
                 # As above, prevent errors if the signal is a single spectrum
674
675
                 output['thickness'] = te
676
677
            return eps, te, Srfint
678
679
680
681
        def im_dielectric_function(self, track_process = False, plot = False):
682
683
             INPUT:
                 self -- the image of which the dielectic functions are calculated
684
685
                 track_process -- boolean, default = False, if True: prints for each pixel
                     that program is busy with that pixel.
                 plot -- boolean, default = False, if True, plots all calculated dielectric
686
                     functions
687
             OUTPUT:
                 self.dielectric_function_im_avg = average dielectric function for each
688
                     pixel
                 self.dielectric_function_im_std = standard deviation of the dielectric
689
                     function at each energy for each pixel
690
                 self.S_s_avg = average surface scattering distribution for each pixel
                 self.S_s_std = standard deviation of the surface scattering distribution at
691
                      each energy for each pixel
692
                 self.thickness_avg = average thickness for each pixel
                 self.thickness_std = standard deviation thickness for each pixel
693
```

```
694
                 self.IEELS_avg = average bulk scattering distribution for each pixel
695
                 self.IEELS_std = standard deviation of the bulk scattering distribution at
                     each energy for each pixel
             .....
696
697
             #TODO
698
             #data = self.data[self.deltaE>0, :,:]
699
             #energies = self.deltaE[self.deltaE>0]
700
             if not hasattr(self, 'ZLPs_gen'):
701
                 self.calc_ZLPs_gen2("iets")
702
             self.dielectric_function_im_avg = (1+1j)*np.zeros(self.data[ :,:,self.deltaE>0
                ].shape)
703
             self.dielectric_function_im_std = (1+1j)*np.zeros(self.data[ :,:,self.deltaE>0
                ].shape)
704
             self.S_s_avg = (1+1j)*np.zeros(self.data[:,:,self.deltaE>0].shape)
705
             self.S_s_std = (1+1j)*np.zeros(self.data[:,:,self.deltaE>0].shape)
706
             self.thickness_avg = np.zeros(self.image_shape)
707
             self.thickness_std = np.zeros(self.image_shape)
708
             self.IEELS_avg = np.zeros(self.data.shape)
709
             self.IEELS_std = np.zeros(self.data.shape)
710
             N_ZLPs_calculated = hasattr(self, 'N_ZLPs')
             #TODO: add N_ZLP saving
711
             #if not N_ZLPs_calculated:
712
                  self.N_ZLPs = np.zeros(self.image_shape)
713
714
            if plot:
715
                 fig1, ax1 = plt.subplots()
716
                 fig2, ax2 = plt.subplots()
717
             for i in range(self.image_shape[0]):
718
                 for j in range(self.image_shape[1]):
719
                     if track_process: print("calculating dielectric function for pixel " ,
720
                     data_ij = self.get_pixel_signal(i,j)#[self.deltaE>0]
721
                     ZLPs = self.calc_ZLPs(i,j)#[:,self.deltaE>0]
722
                     dielectric_functions = (1+1j)* np.zeros(ZLPs[:,self.deltaE>0].shape)
723
                     S_ss = np.zeros(ZLPs[:,self.deltaE>0].shape)
724
                     ts = np.zeros(ZLPs.shape[0])
725
                     IEELSs = np.zeros(ZLPs.shape)
726
                     for k in range(23,28):#ZLPs.shape[0]):
727
                         ZLP_k = ZLPs[k,:]
728
                         N_ZLP = np.sum(ZLP_k)
729
                         IEELS = data_ij-ZLP_k
730
                         IEELS = self.deconvolute(i, j, ZLP_k)
731
                         IEELSs[k,:] = IEELS
732
                         if plot:
733
                             #ax1.plot(self.deltaE, IEELS)
                             plt.figure()
734
                             plt.plot(self.deltaE, IEELS)
735
```

```
736
                         #TODO: FIX ZLP: now becomes very negative!!!!!!
737
                         #TODO: VERY IMPORTANT
738
                         dielectric_functions[k,:], ts[k], S_ss[k] = self.kramers_kronig_hs(
                             IEELS, N_{ZLP} = N_{ZLP}, n = 3)
739
                         if plot:
740
                             #plt.figure()
741
                             plt.plot(self.deltaE[self.deltaE>0], dielectric_functions[k,:]*
                                 2)
742
                             plt.xlim(0,10)
743
                             plt.ylim(-100, 400)
744
                     #print(ts)
745
                     self.dielectric_function_im_avg[i,j,:] = np.average(
746
                         dielectric_functions, axis = 0)
747
                     self.dielectric_function_im_std[i,j,:] = np.std(dielectric_functions,
748
                     self.S_s_avg[i,j,:] = np.average(S_ss, axis = 0)
749
                     self.S_s_std[i,j,:] = np.std(S_ss, axis = 0)
750
                     self.thickness_avg[i,j] = np.average(ts)
751
                     self.thickness_std[i,j] = np.std(ts)
752
                     self.IEELS_avg[i,j,:] = np.average(IEELSs, axis = 0)
                     self.IEELS_std[i,j,:] = np.std(IEELSs, axis = 0)
753
754
             #return dielectric_function_im_avg, dielectric_function_im_std
755
756
        def crossings_im(self):#, delta = 50):
757
             INPUT:
758
759
                 self
             OUTPUT:
760
761
                 self.crossings_E = numpy array (image-shape, N_c), where N_c the maximimun
                     number of crossings of any pixel, 0 indicates no crossing
762
                 self.crossings_n = numpy array (image-shape), number of crossings per pixel
763
             Calculates for each pixel the crossings of the real part of the dielectric
                 function \
764
                 from negative to positive.
765
766
             self.crossings_E = np.zeros((self.image_shape[0], self.image_shape[1],1))
767
             self.crossings_n = np.zeros(self.image_shape)
768
             n_max = 1
769
             for i in range(self.image_shape[0]):
770
                 #print("cross", i)
771
                 for j in range(self.image_shape[1]):
                     #print("cross", i, j)
772
773
                     crossings_E_ij, n = self.crossings(i,j)#, delta)
774
                     if n > n_max:
775
                         #print("cross", i, j, n, n_max, crossings_E.shape)
```

```
776
                          crossings_E_new = np.zeros((self.image_shape[0], self.image_shape[1
777
                         #print("cross", i, j, n, n_max, crossings_E.shape, crossings_E_new
                              [:,:,:n_max].shape)
778
                          crossings_E_new[:,:,:n_max] = self.crossings_E
779
                          self.crossings_E = crossings_E_new
780
                         n_max = n
781
                          del crossings_E_new
782
                     self.crossings_E[i,j,:n] = crossings_E_ij
783
                     self.crossings_n[i,j] = n
784
        def crossings(self, i, j):#, delta = 50):
785
786
             #1 = len(die_fun)
787
             die_fun_avg = np.real(self.dielectric_function_im_avg[ i, j, :])
788
             #die_fun_f = np.zeros(l-2*delta)
789
             #TODO: use smooth?
             11 11 11
790
791
             for i in range(self.l-delta):
792
                 die_fun_avg[i] = np.average(self.dielectric_function_im_avg[i:i+delta])
793
794
             crossing = np.concatenate((np.array([0]),(die_fun_avg[:-1]<0) * (die_fun_avg[1</pre>
                 : ] >= 0)))
795
             deltaE_n = self.deltaE[self.deltaE>0]
796
             \#deltaE_n = deltaE_n[50:-50]
797
             crossing_E = deltaE_n[crossing.astype('bool')]
798
             n = len(crossing_E)
799
             return crossing_E, n
800
801
802
803
        #PLOTTING FUNCTIONS
804
        def plot_sum(self, title = None, xlab = None, ylab = None):
805
806
             INPUT:
807
                 self -- spectral image
808
                 title --- str, delfault = None, title of plot
809
                 xlab -- str, default = None, x-label
                 ylab -- str, default = None, y-label
810
811
             OUTPUT:
             Plots the summation over the intensity for each pixel in a heatmap.
812
813
814
             #TODO: invert colours
815
             if hasattr(self, 'name'):
816
                 name = self.name
817
             else:
                 name = ''
818
```

```
819
             plt.figure()
820
             if title is not None:
821
                 plt.title("intgrated intensity spectrum " + name)
             else:
822
823
                 plt.title(title)
824
             ax = sns.heatmap(np.sum(self.data, axis = 2))
             if not hasattr(self, 'pixelsize'):
825
826
                 plt.xlabel(self.pixelsize)
827
                 plt.ylabel(self.pixelsize)
828
             if xlab is not None:
829
                 plt.xlabel(xlab)
             if ylab is not None:
830
831
                 plt.ylabel = ylab
832
             plt.show()
833
834
835
        def plot_all(self, same_image = True, normalize = False, legend = False,
836
                      range_x = None, range_y = None, range_E = None, signal = "EELS", log =
                           False):
837
             #TODO: add titles and such
838
             if range_x is None:
839
                 range_x = [0, self.image_shape[1]]
840
             if range_y is None:
841
                 range_y = [0,self.image_shape[0]]
842
             if same_image:
843
                 plt.figure()
                 plt.title("Spectrum image " + signal + " spectra")
844
845
                 plt.xlabel("[eV]")
                 if range_E is not None:
846
847
                     plt.xlim(range_E)
848
             for i in range(range_y[0], range_y[1]):
849
                 for j in range(range_x[0], range_x[1]):
850
                     if not same_image:
851
                         plt.figure()
                         plt.title("Spectrum pixel: [" + str(j) +","+ str(i) + "]")
852
853
                         plt.xlabel("[eV]")
854
                         if range_E is not None:
855
                              plt.xlim(range_E)
856
                         if legend:
857
                              plt.legend()
858
                     signal_pixel = self.get_pixel_signal(i,j,signal)
859
                     if normalize:
860
                          signal_pixel /= np.max(np.absolute(signal_pixel))
                     if log:
861
862
                          signal_pixel = np.log(signal_pixel)
863
                         plt.ylabel("log intensity")
```

```
plt.plot(self.deltaE, signal_pixel, label = "[" + str(j) +","+ str(i) +
864
                          "]")
865
                 if legend:
866
                     plt.legend()
867
868
869
        #STATIC METHODS
870
        @staticmethod
871
        def get_prefix(unit, SIunit = None, numeric = True):
             11 11 11
872
             INPUT:
873
874
                 unit -- str, unit of which the prefix is wanted
875
                 SIunit --- str, default = None, the SI unit of the unit of which the prefix
                     is wanted \
876
                              (eg 'eV' for 'keV'), if None, first character of unit is
                                  evaluated as prefix
                 numeric -- bool, default = True, if numeric the prefix is translated to the
877
                      numeric value \
878
                              (e.g. 1E3 for 'k')
879
             OUTPUT:
880
                 prefix -- str or int, the character of the prefix or the numeric value of
                     the prefix
             .....
881
882
             if SIunit is not None:
883
                 lenSI = len(SIunit)
884
                 if unit[-lenSI:] == SIunit:
885
                     prefix = unit[:-lenSI]
886
                     if len(prefix) == 0:
                         if numeric: return 1
887
888
                         else: return prefix
889
                 else:
890
                     print("provided unit not same as target unit: " + unit + ", and " +
                         SIunit)
891
                     if numeric: return 1
892
                     else: return prefix
893
             else:
894
                 prefix = unit[0]
895
             if not numeric:
896
                 return prefix
897
             if prefix == 'p':
898
899
                 return 1E-12
900
             if prefix == 'n':
901
                 return 1E-9
             if prefix == ' ' or prefix == ' ' or prefix == 'u':
902
903
                 return 1E-6
```

```
904
             if prefix == 'm':
905
                 return 1E-3
906
             if prefix == 'k':
907
                 return 1E3
908
             if prefix == 'M':
909
                 return 1E6
             if prefix == 'G':
910
                 return 1E9
911
912
             if prefix == 'T':
913
                 return 1E12
914
             else:
                 print("either no or unknown prefix in unit: " + unit + ", found prefix " +
915
                     prefix + ", asuming no.")
916
             return 1
917
918
919
920
        #CLASS THINGIES
921
        def __getitem__(self, key):
             """ Determines behavior of 'self[key]' """
922
923
             return self.data[key]
924
            #pass
925
926
927
928
        def __str__(self):
             if hasattr(self, 'name'):
929
                 name_str = ", name = " + self.name
930
931
             else:
                 name_str = ""
932
933
             return 'Spectral image: ' + name_str + ", image size:"+ str(self.data.shape[0])
                  + 'x' + \
934
                         str(self.data.shape[1]) + ', deltaE range: [' + str(round(self.
                             deltaE[0],3)) + ',' + \
935
                              str(round(self.deltaE[-1],3)) + '], deltadeltaE: ' + str(round(
                                 self.ddeltaE,3))
936
937
        def __repr__(self):
938
             data_str = "data * np.ones(" + str(self.shape) + ")"
939
             if hasattr(self, 'name'):
                 name_str = ", name = " + self.name
940
941
             else:
942
                 name_str = ""
             return "Spectral_image(" + data_str + ", deltadeltaE=" + str(round(self.
943
                 ddeltaE, 3)) + name_str + ")"
944
```

```
945
         def __len__(self):
946
             return self.l
947
948
949
950
951 def CFT(x, y):
952
         x_0 = np.min(x)
953
         N<sub>0</sub> = np.argmin(np.absolute(x))
954
        N = len(x)
         x_max = np.max(x)
955
956
         delta_x = (x_max - x_0)/N
957
         k = np.linspace(0, N-1, N)
         cont_factor = np.exp(2j*np.pi*N_0*k/N)*delta_x #np.exp(-1j*(x_0)*k*delta_omg)*
958
             delta_x
959
         F_k = cont_factor * np.fft.fft(y)
960
         return F_k
961
962 def iCFT(x, Y_k):
963
         x_0 = np.min(x)
964
         N<sub>0</sub> = np.argmin(np.absolute(x))
965
         x_max = np.max(x)
966
         N = len(x)
967
         delta_x = (x_max - x_0)/N
968
         k = np.linspace(0, N-1, N)
969
         cont_factor = np.exp(-2j*np.pi*N_0*k/N)
970
         f_n = np.fft.ifft(cont_factor*Y_k)/delta_x # 2*np.pi ##np.exp(-2j*np.pi*x_0*k)
         return f_n
971
972
973
974
975
976
977
978
979
980 #%%
981
982
983 #data = np.load("area03-eels-SI-aligned.npy")
984 #energies = np.load("area03-eels-SI-aligned_energy.npy")
985
987
    #dielectric_function_im_avg, dielectric_function_im_std = im_dielectric_function(data,
        energies)
988
```

```
989
 990 #%%
991 #crossings_E, crossings_n = crossings_im(dielectric_function_im_avg, energies)
992
 993
994 #%%
995
 996 #plt.figure()
997 #plt.imshow(crossings_n, cmap='hot', interpolation='nearest')
998 #plt.
999
1000
1001 #ax = sns.heatmap(crossings_n)
1002 #plt.show()
1003
1004 #%%
1005 #dmfile = dm.fileDM('area03-eels-SI-aligned.dm4')
1006 #data2 = dmfile.getDataset(0)
1007
1008 im = Spectral_image.load_data('area03-eels-SI-aligned.dm4')#('pyfiles/area03-eels-SI-
         aligned.dm4')
1009 im.cut_image([0,70], [95,100])
1010 #im.cut_image([40,41],[4,5])
1011 im.calc_ZLPs_gen2(specimen = 4)
1012 im.smooth(window_len=50)
1013 im.im_dielectric_function()
1014 im.crossings_im()
1015
1016 #%%
1017
1018 plt.figure()
1019 plt.title("number of crossings real part dielectric function")
1020 ax = sns.heatmap(im.crossings_n)
1021 plt.show()
1022
1023 plt.figure()
1024 plt.title("energy of first crossings real part dielectric function")
1025 ax = sns.heatmap(im.crossings_E[:,:,0])
1026 plt.show()
1027
1028
1029 plt.figure()
1030 plt.title("thickness of sample")
1031 ax = sns.heatmap(im.thickness_avg)
1032 plt.show()
```